



Physician: J. Laningham, MD  
 Area Treated: Face & Neck  
 Energy Level: 20  
 Sessions: 1  
 Pulse Count: 332

Area Treated	Depth
Cheek	3mm
Neck	3mm
Periorbital	2mm

The treatment technique herein is presented to demonstrate the method utilized by the practitioner. The information on the products and procedures contained in this Case Study does not represent and does not constitute medical advice or recommendations by InMode and should not be relied upon as such. Please note: it is the clinicians' responsibility to decide which application and settings they use in each clinical situation.

